## **RECORD OF RISK ASSESSMENT**

Title of the risk assessment	Equipment Risk Assessment for Metprep Saphir 520 grinding/polishing machine in nCATS Laboratory				
Date risk assessment carried out	30 <sup>th</sup> August 2018				
Describe the work being assessed	Grinding and lapping of samples for surface characterisation and tribological testing for teaching, research and for commercial clients.				

Describe the location at which the work is being carried out	Building 7, room 2023

Where appropriate list the individuals doing the work and the dates/times when the work will be carried out	Visitors, Technical, Academic Staff, Research and Project Students

List any other generic or specific risk	"Risks outside this generic assessment (based on the				
assessments or other documents that	materials employed) will require a separate assessment to be				
relate to this risk assessment – use	made.				
hyperlinks if possible	Undergraduate student (and where appropriate staff)				
	research activities to be risk assessed on individual basis."				

Name and post of risk assessor	Terry Harvey, Area Academic Lead

List the names and post of those assisting in compiling this risk assessment	

Name, post and where required, signature of the responsible manager/supervisor approving the risk assessment	Ling Wang, Head of Group

Reference number and version number	Version One
of risk assessment	

## Assessment

Title of risk assessment

Equipment Risk Assessment for lapping machine in nCATS Laboratory

Γ									Severity						
		Risk Acceptabili	ty		Risk Mat	rix	very Iow	low	medium	high	very high				
	1-3 Risk /	Acceptable					1	2	3	4	5				
	4-6 Risk t	o be reduced if read	ily possible		Certain	5	5	10	15	20	25				ls?
	7-14 Risk t	o be reduced if reas	onably practicable		Likely	4	4	8	12	16	20	p			ontro
	15-25 Risk	Jnacceptable			Possible	3	3	6	9	12	15	lihoo	erity	0	ra co
					Less likely	2	2	4	6	8	10	Overall Likelihood	Overall Severity	Score	or extra controls?
					Improbable	e 1	1	2	3	4	5	erall	erall	lisk 9	es o
												Q	ð	al R	changes
ref	Task/Aspect of wo	k Hazard	Harm and how it could	arise	Who could be affected?		Existing measures to control risk						isk ctors	Residual Risk	Any ch
1	Machinery and equipment	Finger trap; rotating components; sharp edges	Crush injury, bruising, trapping as the sample holder is lowered onto the grinding/lapping discs; cuts bruising and amputation du to entanglement in rotating machinery; cuts on sample that have been ground and lapped producing sharp edges		User				fe operatior			1	3	3	no
2	Fire	Ignition of flammable substances	Burn, smoke inhalatior ignition of flammable chemicals		All	electri	Only small volumes of flammable substances are used, a electrical sources of ignition are away from flammables and insulated				1	2	2	no	
3	Electrical equipment	Electricity	Electrical shock/burn fi contact with mains pov equipment			Installation and maintenance of equipment conducted by qualified electricians. Users will do a quick visual inspection before use. Annual PAT testing.				1	2	2	no		
4	Pneumatic loadin system	g Pressurised air	Injury from pressurised due to escape of compressed air	l jet	User	labora	tories 6 n		andard push pipe. All oth ne rig.			1	2	2	no

## Post Risk Assessment Actions

## Title of risk assessment

Equipment Risk Assessment for lapping machine in nCATS Laboratory

Have any of the specialist control measures listed below been identified as required during risk assessment? – indicate yes or no – if yes then include details on the post assessment action list below.		
Is any exposure monitoring required?	No	
Is any occupational health monitoring required?	No	
Are there any hazards or other factors that could affect pregnant or nursing mothers?	No	

Is any specific training required before people can carry out this work?	Yes
All operators of equipment should have training in that equipment before they carry out any experimental work	

Are there any additional procedures or risk assessments required as a result of this risk assessment?	Yes
Training on test equipment undertaken plus undergraduate students and visitors will be required to complete a Risk Assessment before any testing starts.	

Are there any specialist disposal arrangements required?	No

Are there any special emergency arrangements required?	No

Post Assessment Actions			
Action	By whom	By when	

#### Examples of harm that can Examples of work activities during **Examples of hazards** result if risks are not hazard may be encountered adequately controlled Substances that are harmful if Use or generation during laboratory work, Dermatitis, chemical burn, poisoning contacted, ingested, injected, cleaning activities, outdoor pursuits, or other illness inhaled maintenance work Manual handling lifting, carrying, pushing, pulling, sliding of Bruising, Back injury, strains equipment or people Water watersports, outdoor pursuits, field work, drowning research using flumes Pressure and vacuum systems compressed air or gas systems, vacuum explosion or implosion, injury from rigs pressure jets, hearing damage Psychological working alone, overseas, isolated stress or distress, suicide, long term situations, adverse conditions mental conditions Vehicle moving or manoeuvring vehicles on public Crushing, impact injuries or private roads or yards, towing, cross country Electrical equipment, temporary generators or Electrical shock/burn supplies, experimental rigs, exposed cables, maintenance work Environmental exposure to extremes of heat, cold, wind, Hot burns, cold burns dust during field work or maintenance work Height working at height, outdoor activities Cuts/bruises, Broken bones, Concussion Fire flame cutting equipment, welding or burns, smoke inhalation, brazing, heating equipment, outdoor barbeques or fires Ionising radiation radioactive materials, imaging machines long term illness, burns Machinery and equipment workshop tools, mobile equipment, hand Crushing. trapping, cuts and bruises, tools amputation Non lonising radiation lasers, ultrasound, microwaves surface or deep burns, eyesight damage Noise or vibration agricultural machinery, wind tunnels. hearing loss, hand arm vibration

vehicles. workshop equipment, test rigs

entering tanks, voids in buildings, boilers,

furnaces, sewer and water pipes and

manholes

Confined spaces

Asphyxiation, illness due to breathing harmful gasses or vapours, explosion

syndrome, internal organ damage

### Faculty of Engineering and the Environment

#### Method Statement (Equipment)

Metprep Saphir 520 grinding/lapping machine

Location of Equipment (Building and Room/Laboratory number)	7/2023		Date	30 <sup>th</sup> August 2018
Assessor (Name, ID number)		Contact Details (Email, Telephone number)		
Dr. Terry Harvey, 11467115		harveyt@soton.ac.uk; x23761		
Supervisor		<b>Contact Detail</b> (Email, Teleph	<b>ls</b> one, Room number)	
Prof. Ling Wang		Ling.wang@soton.ac.uk; 7/4081, x25076		

#### Introduction / Overview.

Name of Equipment

(What is the purpose of the equipment? Who is likely to use it?)

The Metprep Saphir 520 grinding/lapping machine is used in the preparation of engineering surfaces for Teaching, Research and Commercial Clients.

A printout of this method statement, the associated Equipment Risk Assessment (ERA), along with a list of users will be kept near the equipment. Also included will be contact details of the equipment leader(s) (person or persons responsible for maintenance and training of the equipment). Also included are a products data sheets for common grinding and lapping consumables, these are classified as not harmful to health or hazardous to the environment and does not present a fire hazard.

Note any activities outside those contained within this document will require an individual risk assessment.

All users will receive a copy of the above documents once they have been trained for unsupervised operation.

#### **Description of Equipment.**

(Provide details of the equipment, what is does and how it does it – the more detail you provide the more likely is anybody reading this will understand what is being done)

The Saphir 520 has a rotating base that the magnet platen is positioned on, see Figure 1. This enables a variety of consumable discs (such MD-Piano, MD-Nap, MD-Allegro, etc.) to be affix to the top by the magnet offering both coarse (grinding) or fine (lapping) removal of material. Lubrication can be applied either by the water spray system. For the cloths (such as MD-Nap) the lapping fluid is applied manually from a dispensing bottle. The amount of lubrication/lapping fluid is varied according the operation and material been processed.




Figure 1 Image of the Saphir 520 grinding/lapping machine with labels for important parts.

During operation the selected grade of grinding disc or lapping cloth disc are place on the magnetic platen, the sample holder, see Figure 2, is lowered into position over the disc and up to five samples are placed in the slots. Using the control panel, see Figure 3 the rotating speed of the base plate (and thus discs) is set. Also the applied load to the sample can be adjusted and grinding/lapping duration can be set. Once the conditions are set, either the lubricant (water) is turned on or lapping fluid applied, then the machine can be started. The grinding/lapping will automatically run for the set period and during this time a load will be applied to each sample via the pneumatic system (internally within the machine) and the base plate (and attached disc) and sample will rotate. During this period the user will have been trained not to approach rotating components and observe operation remotely but close enough to stop the machine if a problem occurs.




Figure 2 Image of the feed system for supplying diamond and lubricant to lapping discs.

Once the rotating components stop and the samples are removed from the holder and either dried if water was used as lubricant or wiped clean, then washed with a solvent (isopropanol) and dried to remove any

trace of the lapping fluid. The samples are then ready for the next step of the grinding/lapping process, which will mean a smaller size of abrasive particle/material; this could either a change of grinding disc or lapping fluid.

Once the grinding/lapping process is finished the machine is cleaned and returned to park position, by raising the sample holder and turning the machine off at the mains socket.



Figure 3 Image of the control panel.

The control panel, see Figure 3, is used to adjusted and select various functions including:

- Adjust the working duration,
- Adjust the load on the samples and whether it is applied individually to the samples or centrally to the sample holder,
- Adjust the rotation speed,
- Adjust the direction of sample holder rotation: options are synchronous (same direction of base unit) or reverse (opposite direction of base unit),
- Select if water lubrication or suspension feed is used, the latter is not available but will be used with the lapping cloths and fluids,
- Save, load and run grinding/lapping programs,
- Turn water on and off,
- Lower and raise sample holder.

During training users are made aware of the health and safety aspects of operation this machine. These include hazards involving the rotating base and sample holder; users are trained to employ the remote operation system, via the control panel, and to remain away from the rig while it is rotating. Also lowering the sample holder into position can lead to finger/hand trapping and users are made aware of this during training.

Other health and safety issues come from the process itself, grinding and lapping can produce sharp edges on samples and users are made aware of this in training. As the machine is electrical powered by single phase mains electricity it is annually PAT tested (and labelled accordingly), users are also trained to perform a visual examination of the mains lead prior to testing. The machine uses air pressure (from main air-line connected to a remote compressor) to pneumatically raise and lower the sample holder, as well as

load the sample holder/samples. The machine has been connected to an existing 6 mm nylon air-line via a push-fit 't'-piece with an appropriate adaptor to the back of the machine. All other pneumatics are contained internal within the machine, for which it is designed.

Identification of risks and risk mitigation

(list all associated risk likely encountered when using the equipment and any existing risk mitigation in place)

Type of Risk: Rotating machinery

Hazards: The machine uses a rotating base and sample holder.

**Mitigation:** The rotation speeds are moderate (100s rpm) and offer no snag points for entanglement. Operation is also remotely performed so that the user should not be near the rotation components during operation. There is an emergency stop button on the front of the machine, which is easily accessible during operation.

Type of Risk: Pressure gas

Hazards: The machine uses compressed air to run the pneumatically powered loading system.

**Mitigation:** The connection to the machine is via standard push-fit employed through the laboratory, all other pneumatic (air pressurised) system are contained within the machine, which is purpose built for this type of operation.

Type of Risk: Electrical equipment

**Hazards:** The lapping machine has an electric operated motor and pump. The motor and its electrical system are contained within the machine. The pump is a self-contained unit, with a lead that connects the main unit.

**Mitigation:** The machine has been designed for purpose and has been supplied by a company. Annual PAT testing ensures that it is safe to use.

Type of Risk: Chemicals

Hazards: Isopropanol solvent is used: is used for general cleaning.

**Mitigation:** Appropriate PPE is supplied and users are required to wear them when handling this chemicals. Waste solvent is collected for hazardous waste disposal. COSHH forms and MSDS are provided.

#### **Control Measures including training, PPE**

(Identify significant hazards and actions/control measures to be taken)

All users will be trained in the safe operation of the polishing machine. Until training is complete all trainees will be supervised and not allowed to operate the machine by themselves.

#### **Emergency Equipment Shutdown Procedure**

(Describe the steps to shut the equipment in the event of an emergency and the location of any emergency stop(s) the equipment has)

In the event of emergency the polishing machine can be stopped using the red emergency stop button on the front of the machine, see Figure 1.

#### Unattended running

(Can this rig be run unattended? And if so what precautions are needed)

Due to the duration of grinding and lapping unattended running is not common and it is recommended that this is not done.

## Faculty of Engineering and the Environment COSHH assessment form

This form must be completed **before** any work with substances hazardous to health is begun, so that a suitable and sufficient assessment of the health risks is made.

Procedure being carried out	Using Saphir 520 grinding/lapping machine			
<b>Location where the substance will be used</b> What supervision or training will the person carrying out the procedure receive?	7/2023 Training by / instruction from facility experienced user	manager of	Review date <sup>1</sup>	30/8/2018
	Name	Signature		Date
Person performing the work	Dr. Terry Harvey			
Supervisor/grant holder	Prof. Ling Wang			
Divisional Safety Officer or other designated person	Dr. Terry Harvey			

<sup>1</sup> This assessment should be reviewed immediately if there is any reason to consider that the original assessment is no longer valid, e.g. due to significant changes in the work activity.

## Attachments

The following documents must be attached:

- Risk assessment identifying the need for the COSHH assessment and clearly indicating the persons potentially at risk (e.g. staff, students, visitors etc.)

- Full description of the procedure.

- MSDS for all substances in 1 a) or b) below
- Any health and safety information provided by supplier in 1 c) below

## 1 Nature of the hazard and risks identified

## a) Chemicals with Health hazards H phrases H300, H301, H304, H310, H311, H314, H318, H330, H331, H334, H340, H341, H350, H351, H360, H361, H370, H371, H372, H373, EUH029, EUH031, EUH032

Name of substance	Hazard phrases	Possible	Risk from single	Risks from	Duration of	Effects could be
	(Refer to MSDS	exposure route	acute exposure	repeated low	adverse effect	hazardous to human
	- must be	(see key		exposure		reproductive systems
	attached)	below) <sup>2</sup>				
Propanol	H319, H336	1, 2	Serious	Not serious	Short term	Not known

 $^{2}(1)$  Contact skin and/or eyes, (2) Inhalation, (3) Injection and/or sharps

# b) Substances with Physical hazards H phrases H200, H201, H202, H203, H204, H205, H220, H221, H222, H223, H224, H225, H226, H228, H240, H241, H242, H250, H251, H252, H260, H261, H270, H271, H272, H280, H281, EUH001, EUH006, EUH014, EUH018, EUH019, EUH044

Name of substance	Hazard phrases (Refer to MSDS - must be attached)	What are the storage requirements for this material? How will they be met?	Quantity used in procedure	Quantity likely to be held in storage	Risk in planned use	Risk in uncontrolled release from storage
Propanol	H225	Supplied bottle or solvent spray bottle	10-100ml	5 litres	Minor	Minor

## c) Substances without a CAS No and no associated H phrases

Name of substance	Nature of the hazard e.g. biological, flammable, explosive, corrosive	Any other information relating to risks arising from this hazard
n/a		

## 2 Use of substance and control of risks

## a) Control measures

Name of substance	Provide a description of the control measures in place to protect the health and safety of both the	List personal
	user and other persons who may be exposed. Control measures should aim to reduce the risks of	protective equipmen
	exposure to the minimum achievable. Consideration should be given to the use of alternative	or containment
	substances which are less hazardous and have a lower risk associated with their use. In this	required
	section you should also provide details of any post reactive products that have been made as a	
	result of the procedure you have followed and the control measures you intend to use to	
	minimise risks associated with these products. Provide details of any monitoring that will be	
	carried out (e.g. for airborne contaminants or of exposed individuals) <sup>3</sup> . (NB: a full description of	
	the procedure must be attached)	
Propanol	Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate	nitrile gloves;
	ventilation. Remove all sources of ignition. Evacuate personnel to safe areas.	laboratory coat;
		safety spectacles

<sup>3</sup> For the majority of work, atmospheric monitoring should not be necessary for protecting health, providing sufficient thought has gone into ensuring the adequacy of control measures in relation to risks, and the control measures are properly used and maintained

## b) Emergency measures

Name of substance	Describe the actions to be taken in the event of uncontrolled release taking into account the quantity of the spill of the substance (i.e. several grams or kilograms), with details of any equipment and/or service required	List equipment and services required
Propanol	Minor spills can be adsorbed with laboratory wipes, very large spills may require evacuation of the laboratory after opening the windows	nitrile gloves; laboratory coat; safety spectacles

## c) Disposal of substance or product resulting from its use.

Name of substance /	Describe the method to be used for disposal of the substance or its products, with details of any	List equipment and
product	control measures, services, labelling, and/or permissions required	services required
Propanol	Disposal initial requires pouring waste liquid into the 'waste solvent bottle', when the bottle is	Hazardous waste
	full it will collected in as 'hazardous waste'	collection



Creation Date 01-Sep-2009

Revision Date 30-May-2018

Revision Number 16

### SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

#### 1.1. Product identification

Product Description: Cat No. :	<u>Propan-2-ol</u> P/7490/08, P/7490/15, P/7490/17, P/7490/21, P/7490/FP21, P/7490/25, P/7490/27, P/7490/DH25, P/7490/MC15, P/7490/PB08, P/7490/PB17, P/7490/PC24, P/7490/PC25, P/7490/21RSS, P/7490/24RSS, P/7490/25RSS, P/7490/34RSS, P/7490/27RSS
Synonyms	2-Propanol; IPA; Isopropyl alcohol; Propan-2-ol; Isopropanol
CAS-No	67-63-0
EC-No.	200-661-7
Molecular Formula	C3 H8 O
Reach Registration Number	01-2119457558-25
1.2. Relevant identified uses of t	he substance or mixture and uses advised against
Recommended Use	Laboratory chemicals.
0	OLIO Industrial was a line of substances as such as is mere particular to destrict sites

Recommended Use	Laboratory chemicals.
Sector of use	SU3 - Industrial uses: Uses of substances as such or in preparations at industrial sites
Product category	PC21 - Laboratory chemicals
Process categories	PROC15 - Use as a laboratory reagent
Environmental release category	ERC6a - Industrial use resulting in manufacture of another substance (use of intermediates)
Uses advised against	No Information available

#### 1.3. Details of the supplier of the safety data sheet

Company E-mail address	Fisher Scientific UK Bishop Meadow Road, Loughborough, Leicestershire LE11 5RG, United Kingdom begel.sdsdesk@thermofisher.com
1.4. Emergency telephone number	Tel: 01509 231166 Chemtrec US: (800) 424-9300 Chemtrec EU: 001 (202) 483-7616

## **SECTION 2: HAZARDS IDENTIFICATION**

#### 2.1. Classification of the substance or mixture

CLP Classification - Regulation (EC) No 1272/2008			
Physical hazards			
Flammable liquids	Category 2 (H225)		
Health hazards			
Serious Eye Damage/Eye Irritation Specific target organ toxicity - (single exposure)	Category 2 (H319) Category 3 (H336)		

#### Propan-2-ol

#### Environmental hazards

Based on available data, the classification criteria are not met

#### 2.2. Label elements



Signal Word

Danger

#### Hazard Statements

H225 - Highly flammable liquid and vapor

H319 - Causes serious eye irritation

H336 - May cause drowsiness or dizziness

#### **Precautionary Statements**

P210 - Keep away from heat/sparks/open flames/hot surfaces. - No smoking

P240 - Ground/bond container and receiving equipment

P261 - Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray

P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

#### 2.3. Other hazards

Substance is not considered persistent, bioaccumulative and toxic (PBT) / very persistent and very bioaccumulative (vPvB)

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

#### 3.1. Substances

Component	CAS-No	EC-No.	Weight %	CLP Classification - Regulation (EC) No 1272/2008
Isopropyl alcohol	67-63-0	200-661-7	>95	Flam. Liq. 2 (H225) Eye Irrit. 2 (H319) STOT SE 3 (H336)

Reach Registration Number	01-2119457558-25

Full text of Hazard Statements: see section 16

#### **SECTION 4: FIRST AID MEASURES**

#### 4.1. Description of first aid measures

Eye Contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.		
Skin Contact	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention if symptoms occur.		
Ingestion	Do not induce vomiting. Obtain medical attention.		
Inhalation	Move to fresh air. Obtain medical attention. If not breathing, give artificial respiration.		
Self-Protection of the First Aider	Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination.		
4.2. Most important symptoms and effects, both acute and delayed			

Breathing difficulties. May cause central nervous system depression: Inhalation of high vapor concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting

#### 4.3. Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically. Symptoms may be delayed.

#### **SECTION 5: FIREFIGHTING MEASURES**

#### 5.1. Extinguishing media

Propan-2-ol

#### Suitable Extinguishing Media

CO<sub>2</sub>, dry chemical, dry sand, alcohol-resistant foam. Cool closed containers exposed to fire with water spray.

#### Extinguishing media which must not be used for safety reasons

Do not use water jet. Do not use a solid water stream as it may scatter and spread fire.

#### 5.2. Special hazards arising from the substance or mixture

Flammable. Risk of ignition. Vapors may form explosive mixtures with air. Vapors may travel to source of ignition and flash back. Containers may explode when heated.

#### **Hazardous Combustion Products**

Carbon monoxide (CO), Carbon dioxide (CO<sub>2</sub>), peroxides.

#### 5.3. Advice for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear. Thermal decomposition can lead to release of irritating gases and vapors.

#### **SECTION 6: ACCIDENTAL RELEASE MEASURES**

#### 6.1. Personal precautions, protective equipment and emergency procedures

Use personal protective equipment. Remove all sources of ignition. Take precautionary measures against static discharges. Avoid contact with skin, eyes and clothing.

#### 6.2. Environmental precautions

Should not be released into the environment. See Section 12 for additional ecological information.

#### Propan-2-ol

#### 6.3. Methods and material for containment and cleaning up

Prevent further leakage or spillage if safe to do so. Remove all sources of ignition. Soak up with inert absorbent material. Take precautionary measures against static discharges. Use spark-proof tools and explosion-proof equipment. Keep in suitable, closed containers for disposal.

#### 6.4. Reference to other sections

Refer to protective measures listed in Sections 8 and 13.

#### **SECTION 7: HANDLING AND STORAGE**

#### 7.1. Precautions for safe handling

Wear personal protective equipment. Keep away from open flames, hot surfaces and sources of ignition. Use explosion-proof equipment. Use only non-sparking tools. Take precautionary measures against static discharges. Do not get in eyes, on skin, or on clothing. Do not breathe vapors or spray mist. To avoid ignition of vapors by static electricity discharge, all metal parts of the equipment must be grounded.

#### Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing before re-use. Wash hands before breaks and at the end of workday.

#### 7.2. Conditions for safe storage, including any incompatibilities

Keep away from heat and sources of ignition. Flammables area. Keep container tightly closed in a dry and well-ventilated place.

#### 7.3. Specific end use(s)

Use in laboratories

### **SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION**

#### 8.1. Control parameters

#### Exposure limits

List source(s): **UK** - EH40/2005 Containing the workplace exposure limits (WELs) for use with the Control of Substances Hazardous to Health Regulations (COSHH) 2002 (as amended). Updated by September 2006 official press release and October 2007 Supplement. **IRE** - 2010 Code of Practice for the Safety, Health and Welfare at Work (Chemical Agents) Regulations 2001. Published by the Health and Safety Authority.

Component	The United Kingdom	European Union	Ireland
Isopropyl alcohol	STEL: 500 ppm 15 min		TWA: 200 ppm 8 hr.
	STEL: 1250 mg/m <sup>3</sup> 15 min		STEL: 400 ppm 15 min
	TWA: 400 ppm 8 hr		Skin
	TWA: 999 mg/m <sup>3</sup> 8 hr		

#### **Biological limit values**

List source(s):

#### Monitoring methods

BS EN 14042:2003 Title Identifier: Workplace atmospheres. Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents.

MDHS70 General methods for sampling airborne gases and vapours

MDHS 88 Volatile organic compounds in air. Laboratory method using diffusive samplers, solvent desorption and gas

#### Propan-2-ol

#### chromatography

MDHS 96 Volatile organic compounds in air - Laboratory method using pumped solid sorbent tubes, solvent desorption and gas chromatography

#### Derived No Effect Level (DNEL) See table for values

Route of exposure	Acute effects (local)	Acute effects (systemic)	Chronic effects (local)	Chronic effects (systemic)
Oral				
Dermal				888 mg/kg
Inhalation				500 mg/m <sup>3</sup>

**Predicted No Effect Concentration** According to our experience and to the information provided to us, the product does not have any harmful effects if it is used and handled as specified. See values below.

Fresh water	140.9 mg/l
Fresh water sediment	552 mg/kg
Marine water	140.9 mg/l
Water Intermittent	140.9 mg/l
Food chain	160 mg/kg
Microorganisms in sewage	2251 mg/l
treatment	
Soil (Agriculture)	28 mg/kg

#### 8.2. Exposure controls

#### **Engineering Measures**

Ensure that eyewash stations and safety showers are close to the workstation location. Use explosion-proof electrical/ventilating/lighting/equipment. Ensure adequate ventilation, especially in confined areas.

Wherever possible, engineering control measures such as the isolation or enclosure of the process, the introduction of process or equipment changes to minimise release or contact, and the use of properly designed ventilation systems, should be adopted to control hazardous materials at source

#### Personal protective equipment

Eye Protection Goggles (European standard - EN 166)

Hand Protection	Protectiv	ve gloves		
Glove material Butyl rubber Nitrile rubber	Breakthrough time > 480 minutes > 360 - 480 minutes	<b>Glove thickness</b> 0.5 mm 0.35 - 0.55 mm	EU standard EN 374	<b>Glove comments</b> Permeation rate < 0.9 μg/cm2/min As tested under EN374-3 Determination of Resistance to Permeation by Chemicals
Viton (R) Neoprene	> 480 minutes < 40 minutes	0.4 mm 0.7 mm		
 Skin and body protection Wear appropriate protective g			ploves and clothing	g to prevent skin exposure

Inspect gloves before use.

Please observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information)

Ensure gloves are suitable for the task: Chemical compatability, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion.

Remove gloves with care avoiding skin contamination.

Respiratory Protection	When workers are facing concentrations above the exposure limit they must use appropriate certified respirators. To protect the wearer, respiratory protective equipment must be the correct fit and be used and maintained properly
Large scale/emergency use	Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced

Revision Date 30-May-2018

Recommended Filter type:Organic gases and vapours filter Type A Brown conforming to<br/>EN14387Small scale/Laboratory useUse a NIOSH/MSHA or European Standard EN 149:2001 approved respirator if exposure<br/>limits are exceeded or if irritation or other symptoms are experienced.<br/>Recommended half mask:- Valve filtering: EN405; Half mask: EN140; plus filter, EN 141<br/>When RPE is used a face piece Fit Test should be conducted

**Environmental exposure controls** No information available.

Propan-2-ol

## **SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES**

#### 9.1. Information on basic physical and chemical properties

Appearance Physical State	Colorless Liquid		
Odor Odor Threshold pH Melting Point/Range Softening Point Boiling Point/Range Flash Point	Alcohol-like No data available 7 -89.5 °C / -129.1 °F No data available 81 - 83 °C / 177.8 - 181.4 °F 12 °C / 53.6 °F	1% aq. sol @ 760 mmHg <b>Method -</b> Abel Closed Cup (BS 2000 Part 170, IP 170, AS/NZS 2106)	
Evaporation Rate Flammability (solid,gas) Explosion Limits	1.7 Not applicable Lower 2 Vol% Upper 12 Vol%	ASTM D 3539 (Butyl acetate = 1.0) Liquid	
Vapor Pressure Vapor Density Specific Gravity / Density Bulk Density Water Solubility Solubility in other solvents Partition Coefficient (n-octanol/wat Component Isopropyl alcohol Autoignition Temperature Decomposition Temperature Viscosity Explosive Properties	43 mmHg @ 20 °C 2.1 @ 20 °C / 68 °F 0.785 Not applicable Miscible No information available	(Air = 1.0) ASTM D-4052 Liquid ASTM E-659 explosive air/vapour mixtures possible Vapors may form explosive mixtures with air	
Oxidizing Properties	No information available		
9.2. Other information			
Molecular Formula Molecular Weight VOC Content(%) Refractive index Surface tension Coefficient of expansion Dielectric constant Heat of vapourisation Specific heat capacity Thermal conductivity	C3 H8 O 60.1 100% (Organic Carbon (by mass) = 59.9 %) (EC/1999/13) 1.377 at 20 °C / 68 °F (ASTM D-1218) 22.7 mN/m at 20 °C / 68 °F 0.0009 / °C 18.6 at 20 °C / 68 °F 665 J/g 3 kJ/kg °C at 20 °C / 68 °F 0.137 W/m °C at 20 °C / 68 °F		

#### SECTION 10: STABILITY AND REACTIVITY

10.1. Reactivity

None known, based on information available

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous PolymerizationHazardous polymerization does not occur.Hazardous ReactionsNone under normal processing.

<u>10.4. Conditions to avoid</u> Heat, flames and sparks. Keep away from open flames, hot surfaces and sources of ignition.

10.5. Incompatible materials

Strong oxidizing agents. Acids. Halogens. Acid anhydrides.

10.6. Hazardous decomposition products

Carbon monoxide (CO). Carbon dioxide (CO<sub>2</sub>). peroxides.

#### **SECTION 11: TOXICOLOGICAL INFORMATION**

#### 11.1. Information on toxicological effects

#### **Product Information**

**FSUP7490** 

(a) acute toxicity; Based on available data, the classification criteria are not met Oral Dermal Based on available data, the classification criteria are not met Inhalation Based on available data, the classification criteria are not met LD50 Oral LD50 Dermal Component LC50 Inhalation 5840 mg/kg (Rat) 13900 mg/kg (Rat) 72.6 mg/L (Rat) 4 h Isopropyl alcohol 12870 mg/kg ( Rabbit ) (b) skin corrosion/irritation; Based on available data, the classification criteria are not met Category 2 (c) serious eye damage/irritation; (d) respiratory or skin sensitization; Respiratory Based on available data, the classification criteria are not met Skin Based on available data, the classification criteria are not met (e) germ cell mutagenicity; Based on available data, the classification criteria are not met (f) carcinogenicity; Based on available data, the classification criteria are not met There are no known carcinogenic chemicals in this product (g) reproductive toxicity; Based on available data, the classification criteria are not met Category 3 (h) STOT-single exposure;

Results / Target organs	Central nervous system (CNS).
(i) STOT-repeated exposure;	Based on available data, the classification criteria are not met
Target Organs	None known.
(j) aspiration hazard;	Based on available data, the classification criteria are not met
Symptoms / effects,both acute and delayed	May cause central nervous system depression: Inhalation of high vapor concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting

## **SECTION 12: ECOLOGICAL INFORMATION**

#### 12.1. Toxicity Ecotoxicity effects

Propan-2-ol

. Do not empty into drains.

Component	Freshwater Fish	Water Flea	Freshwater Algae	Microtox
Isopropyl alcohol	LC50: = 11130 mg/L,	13299 mg/L EC50 = 48	EC50: > 1000 mg/L, 72h	5
	96h static (Pimephales	h	(Desmodesmus	Photobacterium
	promelas)	9714 mg/L EC50 = 24 h	subspicatus)	phosphoreum 5 min
	LC50: > 1400000 µg/L,		EC50: > 1000 mg/L, 96h	
	96h (Lepomis		(Desmodesmus	
	macrochirus)		subspicatus)	
	LC50: = 9640 mg/L, 96h			
	flow-through			
	(Pimephales promelas)			
	,			

## 12.2. Persistence and degradability<br/>PersistenceExpected to be biodegradable<br/>Persistence is unlikely, based on information available.

Bioaccumulation is unlikely

#### 12.3. Bioaccumulative potential

Component	log Pow	Bioconcentration factor (BCF)
Isopropyl alcohol	0.05	No data available

12.4. Mobility in soil Surface tension	The product contains volatile organic compounds (VOC) which will evaporate easily from all surfaces Will likely be mobile in the environment due to its volatility. Disperses rapidly in air 22.7 mN/m at 20 °C / 68 °F
<u>12.5. Results of PBT and vPvB</u> assessment	Substance is not considered persistent, bioaccumulative and toxic (PBT) / very persistent and very bioaccumulative (vPvB).
<u>12.6. Other adverse effects</u> Endocrine Disruptor Information Persistent Organic Pollutant Ozone Depletion Potential	This product does not contain any known or suspected endocrine disruptors This product does not contain any known or suspected substance This product does not contain any known or suspected substance

## **SECTION 13: DISPOSAL CONSIDERATIONS**

#### 13.1. Waste treatment methods

Waste from Residues / UnusedWasteProductson waste

Waste is classified as hazardous. Dispose of in accordance with the European Directives on waste and hazardous waste. Dispose of in accordance with local regulations.

Contaminated Packaging	Dispose of this container to hazardous or special waste collection point. Empty containers retain product residue, (liquid and/or vapor), and can be dangerous. Keep product and empty container away from heat and sources of ignition.
European Waste Catalogue (EWC)	According to the European Waste Catalogue, Waste Codes are not product specific, but application specific.
Other Information	Waste codes should be assigned by the user based on the application for which the product was used. Do not dispose of waste into sewer. Can be incinerated, when in compliance with local regulations.

## **SECTION 14: TRANSPORT INFORMATION**

#### IMDG/IMO

Propan-2-ol

<u>14.1. UN number</u>	UN1219
14.2. UN proper shipping name	Isopropanol (Isopropyl alcohol)
14.3. Transport hazard class(es)	3
14.4. Packing group	Π

ADR

<u>14.1. UN number</u>	UN1219
14.2. UN proper shipping name	Isopropanol (Isopropyl alcohol)
14.3. Transport hazard class(es)	3
14.4. Packing group	II

<u>IATA</u>

14.1. UN number	UN1219
14.2. UN proper shipping name	Isopropanol
14.3. Transport hazard class(es)	3
14.4. Packing group	II

**14.5. Environmental hazards** No hazards identified

14.6. Special precautions for user No special precautions required

14.7. Transport in bulk according to Not applicable, packaged goods Annex II of MARPOL73/78 and the IBC Code

## **SECTION 15: REGULATORY INFORMATION**

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

International Inventories

X = listed.

Component	EINECS	ELINCS	NLP	TSCA	DSL	NDSL	PICCS	ENCS	IECSC	AICS	KECL
Isopropyl alcohol	200-661-7	-		Х	Х	-	Х	Х	Х	Х	Х

#### **National Regulations**

Component	Germany - Water Classification (VwVwS)	Germany - TA-Luft Class
Isopropyl alcohol	WGK 1	

1	Component	France - INRS (Tables of occupational diseases)
	Isopropyl alcohol	Tableaux des maladies professionnelles (TMP) - RG 84

Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment.

#### 15.2. Chemical safety assessment

A Chemical Safety Assessment/Report (CSA/CSR) has been conducted by the manufacturer/importer

#### **SECTION 16: OTHER INFORMATION**

#### Full text of H-Statements referred to under sections 2 and 3

H225 - Highly flammable liquid and vapor

H319 - Causes serious eye irritation

H336 - May cause drowsiness or dizziness

#### Legend

CAS - Chemical Abstracts Service	<b>TSCA</b> - United States Toxic Substances Control Act Section 8(b) Inventory
EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances PICCS - Philippines Inventory of Chemicals and Chemical Substances IECSC - Chinese Inventory of Existing Chemical Substances KECL - Korean Existing and Evaluated Chemical Substances	,
WEL - Workplace Exposure Limit ACGIH - American Conference of Governmental Industrial Hygienists DNEL - Derived No Effect Level RPE - Respiratory Protective Equipment LC50 - Lethal Concentration 50% NOEC - No Observed Effect Concentration PBT - Persistent, Bioaccumulative, Toxic	<ul> <li>TWA - Time Weighted Average</li> <li>IARC - International Agency for Research on Cancer</li> <li>PNEC - Predicted No Effect Concentration</li> <li>LD50 - Lethal Dose 50%</li> <li>EC50 - Effective Concentration 50%</li> <li>POW - Partition coefficient Octanol:Water</li> <li>vPvB - very Persistent, very Bioaccumulative</li> </ul>
<b>ADR</b> - European Agreement Concerning the International Carriage of Dangerous Goods by Road <b>IMO/IMDG</b> - International Maritime Organization/International Maritime Dangerous Goods Code	ICAO/IATA - International Civil Aviation Organization/International Air Transport Association MARPOL - International Convention for the Prevention of Pollution from Ships

Dangerous Goods Code OECD - Organisation for Economic Co-operation and Development BCF - Bioconcentration factor

#### Key literature references and sources for data

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

#### **Training Advice**

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

ATE - Acute Toxicity Estimate

VOC - Volatile Organic Compounds

Use of personal protective equipment, covering appropriate selection, compatibility, breakthrough thresholds, care, maintenance, fit and standards.

First aid for chemical exposure, including the use of eye wash and safety showers.

Fire prevention and fighting, identifying hazards and risks, static electricity, explosive atmospheres posed by vapours and dusts.

Creation Date	01-Sep-2009
Revision Date	30-May-2018
Revision Summary	SDS sections updated, 9.

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

## End of Safety Data Sheet