

SP5 objectives

POSITION	MAGNIFICATION	IMMERSION TYPE	WORKING DISTANCE	NOTES
1	x10	dry	several mm	Phase Contrast, DIC on request**
2	x20	oil, water, glycerol*	240 μm	DIC, UV corrected, correcting collar for immersion medium needs to be set correctly
3	x40	oil	100 μm	DIC on request**
4	x63	oil	100 μm	DIC, UV corrected
5	x100	oil	90 μm	DIC on request**
6	x63	glycerol*	280 μm	DIC, correcting collar for coverslip thickness needs to be set correctly

* 8 parts glycerol, 2 parts dH₂O

** requires removal of stage to fit, notify DAJ in advance